

teii health

Model:U30B

Maximum weight 600lb
(272kg/43 stone)

FCC Warning

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This device and its antenna must not be located or operated in conjunction with any other antenna or transmitter.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure conditions without restriction.

CAUTION

Slippery when wet.



INSTRUCTION MANUAL

Your new electronic personal scale has been designed to accurately measure weight changes (gain or loss) over time and should give you many years of service under normal use.


SPECIFICATIONS

Description	Electronic Weight Scale with voice and 4G transmission (B595-M3)
Display	LCD digital display
Measurement range	Weight 6.6lb - 600lb (3.0kg - 272.0kg) Ambient temperature 14.0°F - 122.0°F (-10°C - 50°C)
Accuracy	Weight ±(0.3+1d) d=1.00g Ambient temperature ±1.8°F (-1°C)
LCD indication	Weight 5 digits display of KG/lb Ambient temperature 3 digits display of F/°C Symbol Low battery/over load indicate (H/L) / Upload status indicate (PAS or FAL)
Power source	1pc 1000mAh rechargeable Lithium battery DC 3.7V or USB power supply (Type-C)
Automatic power off	In 5 seconds
Main unit weight	Approx. 2.56kg (battery included)
Main unit size	L360mm x W360mm x H287mm
LCD size	67.9 x 29mm
Main unit lifetime	10,000 times under normal use
Battery life	Could be used for about 3 months when charge is full (take measurement 3-6 times every day for normal condition)
Accessories	USB cable, instruction manual
Operating environment	Temperature 5~40°C Humidity 15%~93%RH
Storage environment	Air pressure: 98kPa~106kPa Temperature: -20°C~55°C, Humidity: 13%~93%RH Avoid crash, sun burn or rain during transportation

POWER SUPPLY

- Rechargeable lithium battery (1000mAh)

The battery symbol  shows the current battery level, and it will change according to the battery level.

When powered on, the low battery symbol  will be displayed when the battery needs charging along with voice broadcasts "Low battery". And then plug in the device using the included USB cable to charge. Take the USB cable